



MBT 2001P; Molecular Biology and Genetics (Practical)

Molecular Cloning; Basics of Plasmid, Isolation and Purification

Dr Chandresh Sharma

Assistant Professor

Department of Life Science and Biotechnology

Chhatrapati Shahu Ji Maharaj University, Kanpur



Overview



Introduction

Molecular Cloning and Plasmid Vectors

Plasmid Isolation and Purification




Cloning and expression of target gene






Cloning Process

- ✓ Gene of interest is cut out with restriction enzymes (RE)
- ✓ Host plasmid (circular chromosome) is cut with same REs
- ✓ Gene is inserted into plasmid and ligated with ligase
- ✓ New (engineered) plasmid inserted into bacterium (transform)






Cloning (Details)






Cloning (Details)





Protein Expression Bottleneck



- Protein Biochemistry
 - **soluble**, purifiable protein
- Enzymology
 - **soluble**, active protein
 - 0.1-10 mg of protein
- Crystallography
 - **soluble**, crystallizable protein
 - 5-100 mg of protein



Which Vector?

- ✓ Must be compatible with host cell system (prokaryotic vectors for prokaryotic cells, eukaryotic vectors for eukaryotic cells)
- ✓ Needs a good combination of
 - strong promoters
 - ribosome binding sites
 - termination sequences
 - affinity tag or solubilization sequences
 - multi-enzyme restriction site



Plasmids and Vectors

- ✓ Circular pieces of DNA ranging in size from 1000 to 10,000 bases
- ✓ Able to independently replicate and typically code for 1-10 genes
- ✓ Often derived from bacterial “mini” chromosomes (used in bacterial sex)
- ✓ May exist as single copies or dozens of copies (often used to transfer antibiotic resistance)



Key Parts to a Vector

- ✓ Origin of replication (ORI) – DNA sequence for DNA polymerase to replicate the plasmid
- ✓ Selectable marker (Amp or Tet) – a gene, when expressed on plasmid will allow host cells to survive
- ✓ Inducible promoter – Short DNA sequence which enhances expression of adjacent gene
- ✓ Multi-cloning site (MCS) – Short DNA sequence that contains many restriction enzyme sites





Basic elements of a plasmid/vector

pET developed by WF Studier & BA Moffatt in 1986

- 1) **Ap** = ampicillin resistance
- 2) **ori** = ColE1/pBR322 origin of replication
- 3) **lacI** = lac repressor; bind lacO until IPTG induction
- 4) **T7P** = T7 Polymerase promoter
- 5) **lacO** = lac operator where lac repressor binds
- 6) **➡** = multiple cloning site


pET-15b sequence landmarks	
T7 promoter	453-469
T7 transcription start	452
His•Tag coding sequence	362-380
Multiple cloning sites (<i>Nde</i> I - <i>Bam</i> H I)	319-335
T7 terminator	213-259
lacI coding sequence	(866-1945)
pBR322 origin	3882
<i>bla</i> coding sequence	4643-5500



pET-15b cloning/expression region




Gene Introduction (Bacteria)





Bacterial Transformation





Bacterial Transformation


- ✓ Moves the plasmid into bacterial host
- ✓ Essential to making the gene “actively” express the protein inside the cell
- ✓ 2 routes of transformation
 - CaCl_2 + Temperature shock
 - Electroporation
- ✓ Typical transformation rate is 1 in 10,000 cells (not very efficient) for CaCl_2 , but 1 in 100 for electroporation



Electroporator




**25 microfarads = 2500 V
@ 200 ohms for 5 ms**





Electroporation

- Seems to cause disruption in cell membrane
- Reconstitution of membrane leads to large pores which allow DNA molecules to enter
- Works for bacteria, yeast and animal cells







References

- ✓ Protocol: <https://sci-hub.se/10.1101/pdb.prot093344>
- ✓ Youtube: [Lecture 41 : Isolation of Plasmid DNA – YouTube](#)
<https://www.youtube.com/watch?v=Jyk2RzkxUXw>
<https://www.youtube.com/watch?v=04oLyd2mZv8>



Acknowledgement

For Query



chandreshsharma@csjmu.ac.in; sharmac3001@gmail.com